

LABORATORY REPORT

RIQAS®

PATIENT NAME : MRS T H NIROSHA

AGE : 48 Y 0 M 0 D BILL NO : 248980 (Cash [Full])

GENDER : Female REPORT DATE & TIME : Sun 30 Aug, 2026 04:35 PM

SPECIMEN : Serum SAMPLE DATE & TIME : Sun 30 Aug, 2026 03:27 PM

REFERRED BY : Self Requested

CENTER : COOP-HOSPITAL MATARA

TROPONIN I TITRE (QUANTITATIVE)

TEST	RESULT	UNIT	REF. RANGE
------	--------	------	------------

TROPONINE I - TITRE	< 0.02	ng/mL	
---------------------	--------	-------	--

Reference Ranges

0.0 - 0.06ng/ml - Negative

0.06- 0.1 ng/ml - High risk; advisable to repeat the test within 1hr with fresh specimen

Above 0.1 ng/ml - Positive

Comments :

detects lower levels of Troponin I, than conventional Troponin I assays, hence earlier detection of onset.

Troponin I above 0.06 (99th percentile URL) Indicates myocardial necrosis. Troponin I <0.06 does not exclude a clinical diagnosis of acute coronary syndrome (myocardial damage). Suggest repeat quantitative hS Troponin I after 1 to 2 hours to exclude an ongoing myocardial necrosis in highly suspicious cases.

Third Universal Definition of Myocardial infarction, circulation. 2012; 126

Analytical sensitivity : 0.06 ng/mL

Enzyme Immunoassay method specific reference ranges given.

performed using TOSOH AIA - 360 Fully Automated Immunoassay Analyzer

Dr. Thushari Vithanage
MBBS, PG Dip. Path & MD (Chemical Pathology)
Consultant Chemical Pathologist
DGH Matara

Dr. Dileepa Amarasekiriya
MBBS, D. Path, MD (Haematology)
Consultant Haematologist
DGH Matara

Dr. N J Pittigalage
MBBS, Dip. Microbiology, MD. Microbiology
Consultant Microbiologist
DGH Matara

Piyumali

Medical Laboratory Technologist

No 223, Kumarathunga Mawatha, Matara
printed by: 39 (Sun 2026-08-30 05:14pm)

Tel : +94(0)412232893

Manage by MarCoop health care products private Ltd.

Email: cooplabmatara@gmail.com

Web: www.cooperativehospitals.lk